



FDA Registration Number: 3012317718
FDA Owner Operator Number: 10065049

To place an order visit www.fohhealth.com
Questions: email health@foh.cc
For bulk orders email healthorders@foh.cc
or call 305-757-7940, ext.3 for Customer Service

Our Founders




FOH® CEO and Chief Creative Officer, Simone Mayer and Mayda Perez Executive Vice President, manage their business like a family. With over 180 team members world-wide, we are proud to be a Certified Great Place to Work, since 2016 as well as a Certified Minority Owned Business.

Building lasting, collaborative and service driven relationships with our customers is the key to our success and the reason why our customer loyalty and brand reputation is unmatched in the industry.



Company Profile

Company Name	FOH® Inc.
DBAs Brands	Front of the House® FOH Health Essentials™ room360° by FOH®
Mailing Address	7630 Biscayne Boulevard Miami, FL 33138
Telephone Fax	305.757.7940 305.757.7941
Federal Tax ID	75-2991837
FDA Operator No.	10065049
D&B Duns No.	185
Team Member Count	03-647-9023
Type of Business	Manufacturing, Wholesale & Retail NAICS 423990
Incorporated	February 2002



Our Brands

FOH[®]

a global company

Since our start in 2002, we have transformed an industry accustomed to the ordinary, by offering stylishly unexpected and uniquely trend-forward collections for hospitality and foodservice. All Front of the House[®] and room360[°] by FOH[®] solutions are **designed** by us, **manufactured** by us, and **stocked** by us.

In an industry where presentation and customer service are everything, we are proud to have grown organically, through customer word-of-mouth, one chef, GM or hotelier to another, confident that FOH[®] products reflect remarkable craftsmanship, durability, and value.

With the current changes in our Hospitality and Foodservice Industries, FOH[®] is committed to providing our partners with products to ensure they can continue serving their communities safely. With that, we launched FOH Health Essentials[™].



room360[°][®]
by FOH

FOH[™]
Health
ESSENTIALS

FOHTM

Health

ESSENTIALS

As a minority and family owned business, we have also faced hardship from this crisis. As an organization that is known to pivot quickly to meet our customer's needs, we know that offering FOH Health Essentials™ will bring immediate and lasting benefits to our customers and the people they serve and help support our industry.


FOH Health Essentials™ offers a wide range of Personal Protective Equipment (PPE's), including disposable face masks, gowns, gloves, and much more to our partners in need. Our overseas R&D offices, including in-house product development, engineering, and logistics teams, have focused on providing urgent solutions. We currently procure these items, while planning to manufacture our own collections in the future.



Products

Click on page number to go there

Page 7	Non-Contact Infrared Thermometer
Page 8	3-Ply Ear-Loop Face Mask
Page 9	KN95 Respirator
Page 10	N95 Respirator
Page 11	Full Length Protective Face Shield
Page 12	Protective Goggles
Page 13	Disposable Protective Suit
Page 14	45" Disposable Isolation Gown
Page 15	Disposable Nitrile Gloves
Page 16	Disposable HDPE Gloves



Check out our e-commerce site at www.fohhealth.com

Non-Contact Infrared Thermometer

Product Specifications



Product Description

Infrared instant-read digital thermometer.
No-contact forehead reading. For medical use.
Display temperature measurement in °F or °C
within seconds.

Basic Specifications

- Item Number: HTH200WHF22
- Color: White
- Pack Size: 6

Additional Features

- Back-light LCD display
- 2 AA Batteries included
- Product Size: 4" x 1.6" x 6.1" (10.2 x 4.2 x 15.5cm)
- Weight: 3.7oz (107g) Without Battery
- Measurement Range:
 - 86°F (30°C) and 109.4°F (43°C) - Body Temp
 - 32°F (0°C) and 212°F (100°C) - Object Temp
- Accuracy:
 - (89.6°F) 32°C ~ 94.8°F) 34.9°C ± (32.54°F) 0.3°C
 - (95°F) 35°C ~ (107.6°F) 42°C ± (32.36°F) 0.2°C
- Measurement Distance: .78"-1.9" (2-5cm)
- Automatic Shutdown: ≤18 Seconds
- Response Time: 1 second

Quality Standard

Click certification name to view

- FDA listed under product code FLL - Clinical Electronic Thermometer
- CE Compliant under Medical Devices Directive 93/42/EEC
- EN ISO 13485:2016 Certified

Instructions

See pamphlet enclosed with product.



Style may change based on availability

To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

3-Ply Ear-Loop Face Mask



Product Specifications

Product Description

3-layer mask with ear loops. Melt blown cloth and non-woven material.

Basic Specifications

- Item Number: HMA202BLJ2B
- Material: Non-Woven Material
- Color: Blue
- Dimensions: 6.88" x 3.74"
- Pack Size: 50

Additional Features

- Outside Layer: Leak Proof Non-Woven Material 0.88oz (25g)
- Middle Layer: High Density Filter Layer 0.88oz (25g) Melt-Blown Fabric BFE \geq 95%
- Inner Layer: Direct Contact Skin Layer Non-Woven Material
- Length of Earloop: 6 1/2"
- Quad Fold Design, Pleat Depth: 5/8"
- Size: Universal

Quality Standard

Click certification name to view

- FDA listed under product code LYU - Surgical Apparel
- CE Compliant under Medical Devices Directive 93/42/EEC
- ASTM level 1, BFE \geq 95%

Instructions

Not recommended for the Operating Room or other settings using flammable gases. Discard after use.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

KN95 Respirator

Product Specifications



Product Description

KN95 Respirator with ear loops.

Basic Specifications

- Item Number: HMA203WHJ2H
- Material: Melt Blown Cloth and Non-Woven Material
- Color: White
- Pack Size: 20

Additional Features

- Disposable
- Non-Sterile
- 5-Ply
- Outside Layers: 2-Ply Non-Woven Material
- Middle Layer: Melt Blown Cloth
- Inner Layer: Latex Cotton
- Length of Earloop: 8"
- Size: Universal

Quality Standard

Click certification name to view

- FDA listed under product code MSH-Surgical Mask
- CE compliant under European Standard EN149:2001+A1:2009
- ASTM Level 1. BFE \geq 95%. Conforms to KN95, EN149; 2001 FFP2 standards.

Instructions

With sanitized hands, unfold the mask by pulling the sides of the mask apart. Put the mask over your nose and mouth. Adjust the straps to fit comfortably around the ears. Mold the metal bridge around your nose without pinching to make it fit snugly. Perform a seal check before use: cover the mask with both hands without pinching, and inhale and exhale. If you feel air blowing at the edge of the mask, readjust it. Correct fit of the mask requires contact with smooth skin. Discard after use.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

N95 Respirator

Product Specifications



Product Description

4-layer mask with ear loops. Melt blown cloth and non-woven material.

Basic Specifications

- Item Number: HMA215WHJ2H
- Material: Melt Blown Cloth and Non-Woven Material
- Color: White
- Pack Size: 20

Additional Features

- Disposable
- Non-Sterile
- Length of head straps: 10"
- Size: Universal

Quality Standard

Click certification name to view

- FDA listed under product code MSH-Surgical Mask
- CE Compliant under Medical Devices Directive 93/42/EEC
- NIOSH Approved
- BFE \geq 95%
- PFE \geq 95%

Instructions

Position the respirator in your hands with the nose piece at your fingertips. Cup the respirator in your hand allowing the headbands to hang below your hand. Hold the respirator under your chin with the nosepiece up. The top strap goes over and rests at the top back of your head. The bottom strap is positioned around the neck and below the ears. Do not crisscross straps. Place your fingertips from both hands at the top of the metal bridge. Slide fingertips down both sides of the metal strip to mold the nose area to the shape of your nose. To remove: DO NOT TOUCH the front of the respirator! It may be contaminated! Remove by pulling the bottom strap over back of head, followed by the top strap, without touching the respirator. Discard after use.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

Full Length Protective Face Shield

Product Specifications



Product Description

Disposable safety face shield. Full face coverage. Fluid and splatter resistant. Lightweight wrap around design. Elastic band with comfort sponge. Anti-fog and dust proof. Transparent shield.

Quality Standard

Click certification name to view

- FDA listed under product code KPY - Personal Protective Shield
- CE compliant under Directive 2016/425, EN 166:2001

Basic Specifications

- Item Number: HMA200CLT2D
- Full-Length Clear Visor
- Material: Double Anti-Fog PET
- Band Color: Blue
- Dimensions: 13" Width x 8.5" Height
- Shield Thickness: 0.01"
- Sponge Dimension: 10" x 1.5" x 1"
- Pack Size: 200

Instructions

Discard after use. Discontinue use if broken.



Additional Features

- 1" Thick Sponge for Comfort
- Size: Universal

To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

Protective Goggles

Product Specifications



Product Description

Clear PVC frame. High-impact clear anti-fog PC lens. Adjustable strap. Can be worn over eyewear.

Basic Specifications

- Item Number: HMA208CLT2N
- Frame Material: Clear PVC
- Strap Material: Black Elastic
- Dimensions: 7" x 2.25" x 3"
- Pack Size: 15

Additional Features

- Weight: 2.6oz/pc
- Frame Thickness: 0.08"
- Lens Material: High-Impact, Anti-Fog PC
- Lens Thickness: 0.06"
- Size: Universal

Quality Standard

Click certification name to view

- FDA listed under product code HOY-Eye Shield
- CE Mark Certificate under Directive 2016/425, EN 166:2001
- Compliant with ANSI/ISEA Z87.1 – 2015 test standards

Instructions

To clean, dust off any debris and dirt from the goggles, sanitize with soap, and rinse with warm water.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

Disposable Protective Suit

Product Specifications



Product Description

Full body-zip up suit with hood and elastic cuffs. Light weight, coated non-woven material.

Basic Specifications

- Item Number: HGO200WHT2B
- Material: Coated Non-Woven Material
- Color: White
- Pack Size: 50

Additional Features

- Size: Universal
- Chest: 27"
- Sleeve Length: 22"
- Inseam: 30"
- Total Height: 61" (not including hood)
- Non-Sterile. Non-Surgical
- Ankle Style: Elastic
- Cuff Style: Elastic
- Closure Type: Zipper Front with Velcro Panel
- Seam Style: Sewn
- Material Weight: 2.2oz (63g)

Quality Standard

Click certification name to view

- FDA listed under product code OEA -Surgical Apparel
- CE compliant under Medical Devices Directive 93/42/EEC

Instructions

Not recommended for the Operating Room or other settings using flammable gases. Discard after use.



Style may change based on availability

To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

45" Disposable Isolation Gown

Product Specifications



Product Description

45" length gown, Level 1. Light weight and breathable.

Basic Specifications

- Item Number: HGO201YET2M
- Material: SMS/PP/PP+PE+PP Non-Woven Material
- Color: Blue/Green/White/Yellow
- Pack Size: 100

Additional Features

- Length: 45"
- Chest: 54"
- Sleeve length: 22.5"
- Fluid Protection: Minimal
- Non-Sterile. Non-Surgical.
- Size: Universal
- Material Weight: 30g per Cubic Meter

Quality Standard

Click certification name to view

- CE compliant under Medical Devices Directive 93/42/EEC

Instructions

Not Recommended for the Operating Room or other settings using flammable gases. Discard after use.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

Disposable Nitrile Gloves

Product Specifications



Product Description

Latex-free. Powder-free. Box of 100.

Basic Specifications

- Item Numbers:
 - HGL201BLT2F - Small
 - HGL202BLT2F - Medium
 - HGL203BLT2F - Large
 - HGL204BLT2F - X Large
- Material: Nitrile
- Color: Blue
- Pack Size: 10 Boxes of 100

Additional Features

- Sizes: All 9.4" Length (240 mm)
 - Small: 2.75"-3.2" Width (70-81 mm)
 - Medium: 3.2"-3.5" Width (81-90 mm)
 - Large: 3.9"-4.3" Width (100-110 mm)
 - X Large: 4.3"-4.7" Width (110-120 mm)
- Latex Free
- Thickness: .0025" (0.07mm)
- Powder-Free
- Smooth Finish
- Rolled Cuff Style

Quality Standard

Click certification name to view

- FDA listed under product code LZA - Patient Examination Gloves

Instructions

Discard after use.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

Disposable HDPE Gloves

Product Specifications



Product Description

High-density polyethylene gloves. Powder-free. Latex-free. Transparent. One size fits all. Box of 500.

Basic Specifications

- Item Number: HGL230CLT2F
- Material: High-Density Polyethylene
- Color: Clear
- Pack Size: 10 Boxes of 500

Additional Features

- Latex Free
- Powder-Free
- Length: 11.8" (30cm)
- Thickness: .00039" (0.01mm)
- Size: Universal

Quality Standard

Click certification name to view

BRC certified

Instructions

Discard after use.



To place an order visit www.fohhealth.com | For bulk orders email healthorders@foh.cc
Questions: email health@foh.cc or call 305-757-7940, ext.3 for Customer Service

FOH[®]

FOH Health Essentials[™], Front of the House[®] and room360[°] by FOH[®] are divisions of FOH[®], a global healthcare, hospitality and foodservice brand.

